



H3ABioNet

Pan African Bioinformatics Network for H3Africa

16SrRNA Intermediate Bioinformatics Online Course: Int_BT_2019

16S analysis pipeline QC and ASV picking using the dada2 pipeline



H3ABioNet

Pan African Bioinformatics Network for H3Africa



16SrRNA Intermediate Bioinformatics Online Course:
Int_BT_2019
Imane Allali

Outline

- Quality Control
- DADA2 background
- DADA2 workflow

Quality Control

- Before analyzing generated sequence to draw biological conclusions, a quality control check should be performed to make sure there is no biases in the data.
- QC gives a quick impression of whether your data has any problems of which you should be aware before doing any analysis.

Quality Control

Potential problems:

- Low confidence bases (Ns)
- Sequence specific bias
- Sequence contamination
- Adapters
- ...

Quality Control

Software packages for QC:

- FastQC
- MultiQC
- FastX-Toolkit
- PRINSEQ
- TagCleaner
- NGS QC Tool-Kit
- ...

FASTQ format

What is a FastQ file?

FASTQ= FASTA + Quality

FastQ format is a text-based format for storing both a **biological sequence** and its corresponding **quality scores**.



Quality Control

FASTQ format

- Each FastQ file contains hundreds of millions of rows.
- Each **block of 4 lines**, starting with "@" represents a read.

Line 1 begins with a '@' character and is followed by a sequence identifier and an optional description (like a FASTA title line).

Line 2 is the raw sequence letters (ATCG).

Line 3 begins with a '+' character and is optionally followed by the same sequence identifier (and any description).

Line 4 encodes the quality values for the sequence in Line 2, and must contain the same number of symbols as letters in the sequence.

Quality Control

FASTQ format

A FastQ file containing a single sequence might look like this:

@read name

GATTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT

+ read name

!"(((((*+))%%%+)(%%%).1***-+"))**55CCF>>>>>CCCCCCC65

The character '!' represents the **lowest** quality while '~' is the **highest**.

Quality Control

Quality measurements

Base-calling error probabilities are reported by sequencers.
Usually in Phred (quality) score.
Usually coded by ASCII characters

Phred score

$$Q = -10\log_{10}P$$

If the quality of a base is 20, the probability that it is wrong is 0.01

T	C	A	G	T	A	C	T	C	G
40	40	40	40	40	40	40	40	37	35

Quality Control

What is FastQC?

FastQC aims to provide a QC report which can spot problems which originate either in the sequencer or in the library material.

Quality Control

FastQC reports



Normal



Slightly abnormal



Unexpected

FastQC reports


Basic Statistics

Measure	Value
Filename	good_sequence_short.txt
File type	Conventional base calls
Encoding	Illumina 1.5
Total Sequences	250000
Sequences flagged as poor quality	0
Sequence length	40
%GC	45

Quality Control

FastQC reports


Per Base Sequence Quality




FastQC reports

Per Base Sequence Quality

Good quality FastQC report:




Bad quality FastQC report




FastQC reports

Per Sequence Quality Scores

Good quality FastQC report:




Bad quality FastQC report



FastQC reports

Per Base Sequence Content

DNA library



RNA library

